


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	11/16/21	Nebulizing Gas Flow :1.5L/min		:GK11010007	
Time Processed :	7:34:53	CDL Temp :250		Interface Bias :+3.5 kV	
Injection Volume :	0.2	Block Temp :200		Drying Gas Flow :5 L/min	
Sample Name :	Exendin-4			T.Flow :0.2 ml/min	
Cat No. :	RP10874			B.conc :50%H2O/50%MeOH	
Lot No. :	E6963521K				
Theoretical MW :	4186.66				
Observed MW :	4186.4				