


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 10/14/16  
Time Processed : 11:03:04  
Injection Volume : 0.2  
Sample Name : HEL (46-61)  
Cat No. : RP11736  
Lot No. : PE7111610  
Theoretical MW : 1753.83  
Observed MW : 1754.70

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250 °C  
Block Temp : 200 °C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H<sub>2</sub>O/50%MeOH