


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/09/19
Time Processed : 19:00:16
Injection Volume : 0.3
Sample Name : Speract
Cat No. : RP10620
Lot No. : E5564619J
Theoretical MW : 891.90
Observed MW : 891.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH