


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/04/19
Time Processed : 06:13:08 PM
Injection Volume : 0.2
Sample Name : C3a (70-77)
Cat No. : RP20511
Lot No. : PE9031902
Theoretical MW : 823.90
Observed MW : 824.0

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200