


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/25/19
Time Processed : 08:55:13 PM
Injection Volume : 0.3
Sample Name : BMF-BH3
Cat No. : RP20201
Lot No. : E4855019I
Theoretical MW : 3012.50
Observed MW : 3012.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH