


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 03/02/18
Time Processed : 14:48:00
Injection Volume : 0.2
Sample Name : Human Papillomavirus (HPV) E7 protein (49 - 57)
Cat. No. : RP10586
Lot. No. : PE6421801
Theoretical MW : 1120.27
Observed MW : 1120.60

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH