


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 07/25/18
Time Processed : 21:47:46
Injection Volume : 0.2
Sample Name : Simey HY Peptide (738-746)
Cat No. : RP20273
Lot No. : PE2621807
Theoretical MW : 1167.35
Observed MW : 1167.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH