


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/02/17
Time Processed : 23:40:15
Injection Volume : 1
Sample Name : Angiotensin II, human
Cat No. : RP11257
Lot No. : PE8011612
Theoretical MW : 1046.18
Observed MW : 1046.10

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH