


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 02/28/18
Time Processed : 13:05:55
Injection Volume : 0.4
Sample Name : [Sar1, Ile8]-Angiotensin II
Cat No. : RP13720
Lot No. : PE3321802
Theoretical MW : 968.15
Observed MW : 968.70

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH