

Sample Name : β -Amyloid (1-42), human
 Cat. No. : RP 20527
 Lot. No. : PE4881611

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	5	
25.00	Pumps	Pump B Conc.	65	
25.01	Pumps	Pump B Conc.	95	
31.00	Pumps	Pump B Conc.	95	
31.01	Pumps	Pump B Conc.	5	
37.00	Pumps	Pump B Conc.	5	
37.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : UniPS 10-300 10 μ m 4.6 x 250 mm

Chromatogram



1 Det. A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	14.525	79884	10118	0.269
2	14.775	705253	26990	2.378
3	15.710	28357216	931424	95.631
4	17.452	471264	19116	1.589
5	18.378	39114	2426	0.132
Total		29652731	990074	100.000