

Sample Name :TAT peptide
 Cat No. :RP20256
 Lot No. :E7146421K
 Time Processed :16:12:39
 Month-Day-Year Processed :11/17/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	0
5.00	Pumps	B.Conc	0
20.00	Pumps	B.Conc	30
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	0
35.00	Pumps	B.Conc	0
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.859	18384	1625	0.194
2	11.042	9298260	1354113	98.069
3	11.296	143663	17336	1.515
4	11.586	8866	1252	0.094
5	11.989	7249	1165	0.076
6	12.138	4930	841	0.052
Total		9481351	1376333	100.000