


Mass Spectrum



Sample Information

Month-Day Processed : 11/17/21

Time Processed : 11:28:37

Injection Volume : 0.2

Sample Name : TAT peptide

Cat No. : RP20256

Lot No. : E7146421K

Theoretical MW : 1621.91

Observed MW : 1622.4

Interface

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200

Equipment

Interface Bias

Drying Gas Flow

T.Flow

B.conc

: ZJ21010035

: +4.5 kV

: 5 L/min

: 0.2 ml/min

: 50% H₂O / 50% MeOH