


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	05/10/20	Interface Bias	: +4.5kV
Time Processed :	08:59:22 AM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	Thyrotropin-Releasing Hormone (TRH)	B.conc	: 50% H ₂ O / 50% MeOH
Cat No. :	RP10112		
Lot No. :	E5620120E		
Theoretical MW :	362.38		
Observed MW :	362.1		