


Mass Spectrum



Sample Information
Month-Day Processed : 07/03/21
Time Processed : 8:17:34
Injection Volume : 0.4
Sample Name : [Asp371]-Tyrosinase (369-377), human
Cat. No. : RP10270
Lot. No. : E6648621G
Theoretical MW : 1031.17
Observed MW : 1030.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH