


Mass Spectrum



Sample Information
Month-Day Processed : 02/03/21
Time Processed : 4:36:27 PM
Injection Volume : 0.3
Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), human
Cat No. : RP11089
Lot No. : E6222821A
Theoretical MW : 3125.59
Observed MW : 3124.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH