


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/24/17
Time Processed : 17:31:22
Injection Volume : 0.2
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat. No. : RP11110
Lot. No. : PE6771712
Theoretical MW : 2197.60
Observed MW : 2198.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH