


# Mass Spectrum



### Sample Information

User : Gary  
 Month-Day Processed : 06/04/17  
 Time Processed : 13:34:51  
 Injection Volume : 0.2  
 Sample Name : Endothelin 1, human, porcine  
 Cat No. : RPI2758  
 Lot No. : PE4251706  
 Theoretical MW : 2491.90  
 Observed MW : 2492.40

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250°C  
 Block Temp : 200°C

Interface Bias : +4.5kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2ml/min  
 B.conc : 50% H2O/50% MeOH