


Mass Spectrum



Sample Information

Acquired by: Gary
Month-Day Processed: 04/13/20
Time Processed: 09:26:53 AM
Injection Volume: 0.3
Sample Name: Brain Natriuretic Peptide (BNP) (1-32), human
Sample ID: U322JFC250-1
Theoretical MW: 3464.10
Observed MW: 3463.5

Interface: ESI
Nebulizing Gas Flow: 1.5L/min
CDL Temp: 250
Block Temp: 200

Equipment: GK11010007
Interface Bias: +4.5 kV
Drying Gas Flow: 5 L/min
T.Flow: 0.2 ml/min
B.conc: 50% H2O/50% MeOH