


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/17/19
Time Processed : 10:01:43 PM
Injection Volume : 0.3
Sample Name : IRBP(161-180)
Cat No. : RP11841
Lot No. : E48793191
Theoretical MW : 2209.51
Observed MW : 2208.9

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200