


Mass Spectrum



Sample Information
Month-Day Processed : 01/04/21
Time Processed : 7:35:29 PM
Injection Volume : 0.3
Sample Name : Brain Natriuretic Peptide (BNP) (1-32), human
Cat No.: RP11119
Lot No.: E7838220L
Theoretical MW : 3464.10
Observed MW : 3463.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH