


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/27/17
Time Processed : 11:52:57
Injection Volume : 0.2
Sample Name : Neurotensin
Cat. No. : RP10439
Lot. No. : PE6811710
Theoretical MW : 1671.92
Observed MW : 1672.60

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH