


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 04/13/20
Time Processed : 11:00:46 PM
Injection Volume : 0.3
Sample Name : Somatostatin 28 (1-14)
Cat No. : RP10230
Lot No. : E3955320D
Theoretical MW : 1528.70
Observed MW : 1528.5

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200