


Mass Spectrum



Sample Information
Month-Day Processed : 04/21/21
Time Processed : 22:56:02
Injection Volume : 0.4
Sample Name : Jagged - 1 (188 - 204)
Cat No. : RP20331
Lot No. : E9396821D
Theoretical MW : 2107.40
Observed MW : 2106.9

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH