

Sample Name :Endothelin 1, human, porcine
 Cat No. :RP12758
 Lot No. :E3603321E
 Time Processed :10:45:34 AM
 Month-Day-Year Processed :05/18/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>


| Time | Module | Command | Value |
|-------|------------|---------|-------|
| 0.01 | Pumps | B.Conc | 5 |
| 25.00 | Pumps | B.Conc | 65 |
| 25.01 | Pumps | B.Conc | 95 |
| 32.00 | Pumps | B.Conc | 95 |
| 32.01 | Pumps | B.Conc | 5 |
| 40.00 | Pumps | B.Conc | 5 |
| 40.01 | Controller | Stop | |

<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ19010324

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

| Peak# | Ret. Time | Area | Height | Area% |
|-------|-----------|---------|---------|---------|
| 1 | 14.312 | 31059 | 5734 | 0.416 |
| 2 | 14.619 | 28540 | 5100 | 0.382 |
| 3 | 15.770 | 85025 | 16442 | 1.139 |
| 4 | 16.177 | 106978 | 10431 | 1.433 |
| 5 | 16.555 | 7096719 | 938072 | 95.074 |
| 6 | 16.888 | 49080 | 14016 | 0.658 |
| 7 | 18.529 | 28145 | 4925 | 0.377 |
| 8 | 24.272 | 38841 | 5432 | 0.520 |
| Total | | 7464388 | 1000153 | 100.000 |