


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 05/07/17
 Time Processed : 8:28:48
 Injection Volume : 0.2
 Sample Name : Calcitonin Gene Related Peptide (CGRP) II, rat
 Cat No. : RP11092
 Lot No. : PE7511704
 Theoretical MW : 3805.40
 Observed MW : 3804.40

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250 °C
 Block Temp : 200 °C

Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH