


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 01/12/18  
Time Processed : 16:00:08  
Injection Volume : 0.2  
Sample Name : Cys-TAT(47-57)  
Cat. No. : RP20343  
Lot. No. : PE7661801  
Theoretical MW : 1661.99  
Observed MW : 1662.00

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias  
Drying Gas Flow  
T.Flow  
B.conc

: +4.5 kV  
: 5 L/min  
: 0.2 ml/min  
: 50%H<sub>2</sub>O/50%MeOH