


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 10/24/18
Time Processed : 17:43:17
Injection Volume : 0.2
Sample Name : [Sar¹, Ile⁸-⁷-Angiotensin II
Cat. No. : RP13720
Lot. No. : PE4541810
Theoretical MW : 968.15
Observed MW : 968.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH