


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/10/19
Time Processed : 05:06:01 PM
Injection Volume : 0.4
Sample Name : Enfuvirtide (T-20)
Cat No. : RP17307
Lot No. : E4907119I
Theoretical MW : 4491.88
Observed MW : 4491.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH