


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 07/14/17
Time Processed : 11:14:58
Injection Volume : 0.2
Sample Name : Gastrin-1, human
Cat No. : RP12740
Lot No. : PE3301707
Theoretical MW : 2098.21
Observed MW : 2097.70

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH