


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 02/14/18
Time Processed : 12:46:49
Injection Volume : 0.2
Sample Name : TAT
Cat No. : RP12784
Lot No. : PE5731802
Theoretical MW : 1558.85
Observed MW : 1558.95

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH