


Mass Spectrum



Sample Information
Month-Day Processed : 12/26/20
Time Processed : 7:41:43 PM
Injection Volume : 0.4
Sample Name : Protease-Activated Receptor-2, amide
Cat No. : RP20334
Lot No. : E7212820L
Theoretical MW : 614.80
Observed MW : 614.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH