


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	02/01/20	Interface Bias	: +4.5kV
Time Processed :	07:27:55 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	Bradykinin (1-7)	B.conc	: 50%H ₂ O/50%MeOH
Cat No. :	RP11133		
Lot No. :	E7238320B		
Theoretical MW :	756.90		
Observed MW :	756.8		