


Mass Spectrum



Sample Information	Equipment	: GK11010007
Month-Day Processed :	Interface Bias	: +4.5 kV
Time Processed :	Drying Gas Flow	:5 L/min
Injection Volume :	T.Flow	:0.2 ml/min
Sample Name :	B.conc	:50%H2O/50%MeOH
Cat No. :		
Lot No. :		
Theoretical MW :		
Observed MW :		

Interface	: ESI
Nebulizing Gas Flow	: 1.5L/min
CDL Temp	: 250
Block Temp	: 200