


# Mass Spectrum



Sample Information  
User : Gary  
Month-Day Processed : 09/18/19  
Time Processed : 08:08:46 PM  
Injection Volume : 0.3  
Sample Name : Autocamtide 2-amide  
Cat No. : RP20277  
Lot No. : E4872019I  
Theoretical MW : 1526.79  
Observed MW : 1526.4

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50% $\text{H}_2\text{O}$ /50% $\text{MeOH}$