


Mass Spectrum



Sample Information

Month-Day Processed :	11/13/21	Equipment	: GK11010007
Time Processed :	20:17:19	Interface Bias	: +4.5 kV
Injection Volume :	0.5	Drying Gas Flow	: 5 L/min
Sample Name :	β -Endorphin, human	T.Flow	: 0.2 ml/min
Cat No. :	RP11344	B.conc	: 50% H ₂ O/50% MeOH
Lot No. :	E7126921K		
Theoretical MW :	3465.10		
Observed MW :	3464.5		