

Sample Name :Endothelin 1, human, porcine
 Cat. No. : RP12758
 Lot. No. :E3835149H
 Time Processed :13:38:48
 Month-Day-Year Processed :08/16/2019

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	5	
25.00	Pumps	Pump B Conc.	65	
25.01	Pumps	Pump B Conc.	95	
27.00	Pumps	Pump B Conc.	95	
27.01	Pumps	Pump B Conc.	5	
35.00	Pumps	Pump B Conc.	5	
35.01	Controller	Stop		


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment :SS-CM-0310

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	13.767	1187	106	0.031
2	17.445	5445	579	0.140
3	17.673	4353	557	0.112
4	17.947	60639	7354	1.561
5	18.226	50465	4661	1.299
6	18.493	3741832	590762	96.345
7	18.758	10155	1535	0.261
8	19.093	7170	1324	0.185
9	19.526	1452	196	0.037
10	23.162	1070	154	0.028
Total		3883769	607228	100.000