


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 11/25/17
 Time Processed : 11:31:13
 Injection Volume : 0.2
 Sample Name : Pep-1 (uncapped)
 SCat. No. : RP10130
 Lot. No. : PE5291711
 Theoretical MW : 2848.30
 Observed MW : 2848.60

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH