


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/13/17
Time Processed : 22:51:39
Injection Volume : 1
Sample Name : SLLK, Control Peptide for TSP1 Inhibitor
Cat No. : RPI9989
Lot No. : PE8221701
Theoretical MW : 459.59
Observed MW : 459.50

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH