


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/30/20
Time Processed : 10:43:23 AM
Injection Volume : 0.3
Sample Name : HA Peptide
Cat No. : RP11735
Lot No. : E1986120E
Theoretical MW : 1102.15
Observed MW : 1101.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH