


Mass Spectrum



Sample Information

| | | | |
|-----------------------|-------------|-----------------|------------------|
| User : | Gary | Equipment | : GK11010007 |
| Month-Day Processed : | 09/18/19 | Interface Bias | : +4.5kV |
| Time Processed : | 03:27:02 PM | Drying Gas Flow | : 5 L/min |
| Injection Volume : | 0.3 | T.Flow | : 0.2ml/min |
| Sample Name : | C3a (70-77) | B.conc | : 50%H2O/50%MeOH |
| Cat No. : | RP20511 | | |
| Lot No. : | E4942019I | | |
| Theoretical MW : | 823.90 | | |
| Observed MW : | 823.8 | | |