


Mass Spectrum



Sample Information
 Month-Day Processed : 11/10/20
 Time Processed : 09:15:01 AM
 Injection Volume : 0.3
 Sample Name : Gastric Inhibitory Peptide (GIP), human
 Cat No. : RP10795
 Lot No. : E6152520J
 Theoretical MW : 4983.60
 Observed MW : 4983.6

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH