


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/02/19
Time Processed : 05:10:55 PM
Injection Volume : 0.2
Sample Name : HA Peptide
Cat No. : RPI1735
Lot No. : PE8431902
Theoretical MW : 1102.15
Observed MW : 1101.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH