


Mass Spectrum



Sample Information
 Month-Day Processed : 10/23/20
 Time Processed : 20:56:12
 Injection Volume : 0.3
 Sample Name : Neuropeptide Y(1-36)
 Cat No : RP20520
 Lot No : E2269520J
 Theoretical MW : 4272.70
 Observed MW : 4272.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H2O/50% MeOH