


Mass Spectrum



Sample Information
Month-Day Processed : 04/30/21
Time Processed : 7:28:32 AM
Injection Volume : 0.7
Sample Name : IKKY NBD Inhibitory Peptide
Cat. No. : RP20478
Lot. No. : E2007221D
Theoretical MW : 3693.30
Observed MW : 3692.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH