


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 04/24/18
Time Processed : 12:05:46
Injection Volume : 0.2
Sample Name : Ghrelin_rat
Cat. No. : RP10781
Lot. No. : PE8431804
Theoretical MW : 3314.80
Observed MW : 3315.25

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH