


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 12/09/16
Time Processed : 10:50:48
Injection Volume : 0.2
Sample Name : HIV-1 TAT (48-60)
Cat No. : RP20423
Lot No. : PE4681612
Theoretical MW : 1719.00
Observed MW : 1719.20

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH