


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/18/18
Time Processed : 16:03:54
Injection Volume : 0.2
Sample Name : C-Peptide (3-33), human
Cat No. : RPI1116
Lot No. : PE8021810
Theoretical MW : 3020.28
Observed MW : 3021.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH