


Mass Spectrum



Sample Information

User : Gary

Month-Day Processed : 09/17/19

Time Processed : 22:13:15

Injection Volume : 0.3

Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), human

Cat No. : RP11089

Lot No. : E4795219I

Theoretical MW : 3125.59

Observed MW : 3125.5

Equipment : GK11010007

Interface Bias : +4.5kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2ml/min

B.conc : 50%H2O/50%MeOH

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200