


Mass Spectrum



Sample Information

Month-Day Processed :	06/10/21	Equipment	: GK11010007
Time Processed :	4:39:52 PM	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	: 5 L/min
Sample Name :	PLP(139-151)	T.Flow	: 0.2 ml/min
Cat. No. :	RP10550	B.conc	: 50%H2O/50%MeOH
Lot. No. :	E1322221F		
Theoretical MW :	1521.74		
Observed MW :	1521.3		