


Mass Spectrum



Sample Information
Month-Day Processed : 10/15/20
Time Processed : 6:24:03 PM
Injection Volume : 0.4
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat. No. : RP11110
Lot. No. : E2474720J
Theoretical MW : 2197.60
Observed MW : 2197.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH